

In a class by itself, the SWAT-T™ is a multi-function trauma care device.

**Innovation and Versatility in  
TRAUMA CARE**

# SWAT-T™

## STRETCH-WRAP-AND-TUCK TOURNIQUET



### RESEARCH

**2013:** Tourniquets and occlusion: the pressure of design.

Military Medicine : May 2013

*Conclusion: The SWAT-T™ performed better than the CAT (Combat Application Tourniquet) with safer pressures, and was more effective.*

**2012:** Lighting did not affect self-application of a stretch and wrap style tourniquet.

Journal of Special Operations Medicine: Fall 2012

*Conclusion: The SWAT-T™ stretch and wrap style tourniquet can be self-applied properly even in darkness. When properly applied, it stops limb arterial flow.*

**2012:** Stretch and wrap style tourniquet effectiveness with minimal training.

Military Medicine: November 2012

*Conclusion: Proper application of the SWAT-T™ is easy and will stop extremity arterial flow but requires some training for many applicers.*

**2010:** An Evaluation of the SWAT-T

Presented at the 2010 Annual meeting of the Special Operations Medical Association - Douglas M. Kleineir, PhD

*Conclusion: The SWAT-T™ rated very well in the variables evaluated and was rated comparable to, or better than the C-A-T for size, versatility intuitiveness, durability cost and comfort.*



In a class by itself, the SWAT-T™ is not simply a tourniquet but is a trauma multi-tool - with multiple documented saves. Designed by a former Operator/Medic, the SWAT-T will treat a variety of life-threatening injuries. The SWAT-T™ is being carried by EMS, Law Enforcement Officers, as well as Contractors, Federal and Military (Conventional and Special Operations medical and non-medical personnel).

**“The KISS Principle on Steroids!”**

Rob Prevost  
Operations Manager - Triple Canopy

**“100% effective at stopping arterial blood flow... I give it 4.99 out of 5 stars... in the tactical environment, nothing is perfect.”**

Dr. Sydney Vail  
Trauma Surgeon Instructor, International School of Tactical Medicine  
(after independent testing in his vascular lab)

**“Wider tourniquet is more effective for large thigh application and results in less pain... small profile, multifunction, and durability... a great combination!”**

Dr. Matthew Sena  
Trauma Surgeon, Former Special Forces 18D

TRAINING VIDEOS AT  
**SWAT-T.com**



GSA Available  
NSN Pending  
Patents Pending

### SWAT -Tourniquet™

Multi-function tourniquet: Pressure dressing & All-purpose wrap  
Lower space, weight, and cost  
Intuitive design = less training needed!

