

MedEvac4 STANAG 2040 NATO Stretcher



The MedEvac4 Fold STANAG 2040 Stretcher is the NATO version of the reliable and proven MedEvac4 Stretcher Platform, and is manufactured in accordance with NATO Standards (STANAG 2040). It folds into 4 for compact transport / storage, and yet maintains the extremely high load bearing capability that our stretchers are renown for.

CE This appliance conforms to the Directive 93/42/CEE “Medical Devices”

| Specifications: | |
|----------------------|---|
| Dimensions Extended | 2290mm X 584mm X 171mm |
| Dimensions Folded | 630mm X 215mm X 220mm |
| Canvas Specification | Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm Compliant with MIL-L-16446A Standards |
| Frame | Aerospace Grade Aluminum |
| Handles | Retractable, Automatic locking, Ergonomically Designed |
| Weight | 7kg |
| CE Mark | Yes |

*Specifications are subject to change without prior notice.